

# HANS RUDOLPH, INC. 7600 SERIES V2 MASK™ ORO-NASAL CPAP/BILEVEL REUSABLE with VENTS & AAV INSTRUCTIONS FOR USE


**Figure 1**

**Figure 2**


**Figure 3**


**Figure 4**

**Figure 5**

**Figure 6**

## Intended Use:

The 7600 V2 Mask is a reusable, multi-patient, multi-use, Oro-Nasal CPAP/BiLevel mask which incorporates a passive, continuous flow exhaust port. It is intended for use on adults (>30 kg) in homes, hospitals and other clinical settings with CPAP/BiLevel machines for obstructive sleep apnea treatment, and for use with other similar ventilators that use this mask exhaust port configuration providing minimum 3 cm H2O pressure at mask.

## Caution:

Federal law restricts this device to sale by or on physician's order.

## Warnings:

1. Shipped clean, nonsterile. If disinfected or sterile use is required, follow the procedures in this document prior to use of the mask.
2. To avoid skin irritation, never use any cleaning solutions or chemicals on this mask or headgear other than a mild detergent and water.
3. Patients with facial hair may experience mask leakage.
4. Do not block mask vent holes or patient circuit flow.
5. Do not use with ventilators that require double-limb patient circuits.
6. Not for use with CPAP/BiLevel devices which include separate patient circuit exhaust port.
7. When CPAP/BiLevel device is turned off, or is not operating properly, flow through mask vent holes may be inadequate to clear all exhaled gas and some rebreathing may occur which could lead to suffocation.
8. Discontinue mask use if irritation or allergic reaction develops due to mask and contact a healthcare professional.
9. Patient failure to arouse and remove mask after vomiting could result in vomitus aspiration.
10. If supplemental oxygen is used, inhaled percent oxygen will vary depending upon pressure settings, patient breathing pattern, mask size & leak rate.
11. Oxygen flow must be turned off when CPAP/BiLevel system is not operating.
12. Combination of this mask with devices other than for which mask is intended can alter mask safety and performance becoming unsafe or life-threatening.
13. Do not clean using an automated washer-disinfector as this process may remove device markings or finishes.
14. Frequency of cleaning, methods of cleaning or the use of cleaning agents, other than those specified in this document, or exceeding the number of processing cycles can have an adverse effect on the mask and consequently the safety or the quality of the therapy.
15. Mask fall into the semi-critical device category based on potential risk of infection. Although sterilization is recommended whenever practical, high level disinfection should be performed between uses on different patients.

## Contraindications:

1. Open wounds prone to infection.
2. Hemodynamic or cardiorespiratory instability.
3. Unconsciousness or patients unable to remove mask.
4. Claustrophobia, anxiety, or other discomfort with Oro-Nasal mask.
5. Facial or nasopharyngeal deformity or other inability to fit mask & seal properly.
6. Excessive reflux, GI blood or other secretions.
7. Impaired cough reflex, hiatal hernia, or inability to swallow or clear secretions.
8. Barotrauma, upper airway obstruction or facial trauma.

9. Need for ventilation or ventilator support > 12 hours per day.
10. Recent facial, esophageal, or gastric surgery.
11. Patients under medication with drug that may cause vomiting.
12. Patients requiring immediate intubation.

## Risks and Complications:

1. Biosensitivity to device materials (irritation, sensitivity, or other reaction),
2. Irritation/abrasion/ulceration/pain/deformity/drying (face, head, neck, eyes, mucosa, dental),
3. Inhaled, aspirated, applied or handled contaminants (particulates, microbes, body fluids, chemicals),
4. Infection (skin, eyes, mucosa, respiratory),
5. Excessive rebreathed dead space or breathing resistance (dyspnea, atelectasis, barotrauma),
6. Excessive/damaging sound level,
7. Low secretion clearance (especially during upper respiratory infections),
8. Inadequate or excessive mask operating pressure (barotrauma, dyspnea, atelectasis, gastric distention, flatulence),
9. Adverse interactions with attached non-HRI accessories.

## Safety Information:

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

| SYMBOLS                                        |                                        |
|------------------------------------------------|----------------------------------------|
| Medical device manufacturer                    | Needs protection from UV light sources |
| Do not use if the package is damaged or opened | Temperature limitations                |
| Use-by date (first use)                        | Consult the instructions for use       |
| Manufacturer's batch/lot code                  | European conformity (CE) mark          |
| Manufacturer's catalogue/part number           | Indicates product is a medical device  |
| Date of the manufacture                        | Model (Series) number                  |
| For sale by/on physician's order only (USA)    | Size large                             |
| Size medium                                    | Size small                             |
| Size extra small                               | Size petite                            |
| United Kingdom conformity mark                 |                                        |

## Mask Components (fig 1):

1. Mask face piece has plastic braces with slot openings for attachment of headgear strap mounting quick-release clips.
2. Swivel port assembly comprises of sampling port, micro vent holes, anti-asphyxia valve (AAV) and detachable 22 mm (ISO 5356-1) swivel port.
3. AAV (Fig. 2) snaps into swivel port assembly elbow using "locking tabs" for engagement.

## Application & Fitting (fig 4):

1. Determine patient mask size with optional Mask Sizing Gauge (fig 3).
2. Place mask over patient's nose and mouth (1). Chin should fit into mask chin cup portion with top of mask sealing area on the bridge of the nose.
3. Slide headgear over patient's head (2) with one of the lower quick release strap clips disconnected. Reconnect headgear strap clip after headgear is placed (3). Bottom straps should be positioned below ears and top straps above ears & below eyes.
4. Adjust headgear Velcro™ strap tension by pulling back slightly on both the top and bottom straps (4 & 5). Adjust straps at crown of head by pulling each in opposing directions until reaching required tension and then reconnecting Velcro™ (6).
5. Connect mask to patient circuit by following ventilator operating instructions.
6. Turn on ventilator. If you detect leaks reposition mask or adjust strap tension. If leak continues, regardless of adjustments, try another mask size.

Safety Features Verification:

- 1. AAV (fig 5) allows room air breathing for patient if ventilator device stops. AAV diaphragm CLOSES (>3 cm H2O) the room air port in mask during CPAP/ BiLevel ventilation. It OPENS (<3 cm H2O) mask to room air breathing when ventilator device is off. Mask should not be used if AAV is missing, damaged or not functioning properly.
- 2. Mask Removal with Quick-Release Headgear (fig 6). To remove mask slide one of the lower strap quick-release clips slightly forward towards mask front. It will disconnect from mask; pull mask and headgear to opposite side or up over patient's head.

Disassembly for Cleaning:

- 1. Separate headgear from mask.
- 2. Separate swivel port assembly from mask, AAV from swivel port (squeeze two locking tabs & pull AAV out) and rubber cap plug from mask adapter sampling port.

Cleaning Mask face piece, swivel port components & headgear:

Manual Cleaning (Mask & swivel port assembly):

- 1. Soak components for 5 minutes in warm water with mild detergent (neutral pH).
- 2. Hand wash components with sponge or soft brush; then rinse with warm tap water.
- 3. Place components in warm water bath and agitate for two minutes.
- 4. Rinse in clean tap water for at least 1 minute.
- 5. Allow components to air dry or dry with clean, lint free cloth.

Manual Cleaning (Headgear):

- 1. Submerge headgear in warm soapy tap water; gently rub all areas.
- 2. Rinse in clean tap water for 1 minute or until all soap is removed; air dry.
- 3. Precaution: Do not use bleach, chlorine or alcohol based solutions to clean mask or headgear. These solutions can damage components. Direct sunlight exposure on mask and headgear components can cause deterioration and reduce product life.

High Level Disinfection (Mask & Swivel port components less Headgear)

|                                                                                                                                                                                                        |                                                 |                                                                                      |
|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------|--------------------------------------------------------------------------------------|
| Cidex Liquid Glutaraldehyde                                                                                                                                                                            | Steris Revital-Ox Resert                        | Hot Water Pasteurization                                                             |
| Thoroughly clean mask prior to disinfection or sterilization                                                                                                                                           |                                                 |                                                                                      |
| Soak mask components as per Cidex instructions                                                                                                                                                         | Soak mask components as per Steris instructions | Immerse mask components in hot water temperatures of 71.1°C to 76.6°C for 30 minutes |
| Rinse 1: Soak mask components in 1500 ml of sterile water for ≥ 1 minute                                                                                                                               |                                                 |                                                                                      |
| Rinse 2: Soak mask components in 1500 ml of sterile water for ≥ 1 minute                                                                                                                               |                                                 |                                                                                      |
| Dry mask components with clean lint free cloth (preferably sterile)                                                                                                                                    |                                                 |                                                                                      |
| Visually inspect all components for cleanliness, function and defects. If there are any residues, stains or organic debris then repeat the previous steps. Dispose of and replace all defective parts. |                                                 |                                                                                      |

Sterilization (Face Mask assembly less Headgear)

- 1. Thoroughly clean mask prior to disinfection or sterilization
- 2. Mask should be completely assembled (less headgear).

Steam Sterilization

| Parameter          | Steam (Pre-Vacuum)            | Steam (Gravity Displacement) |
|--------------------|-------------------------------|------------------------------|
| Temperature        | 132.2 +3/-1°C                 | 132.2 +3/-1°C                |
| Sterilization Time | 4 minutes                     | 30 minutes                   |
| Dry Time           | 10 minutes                    | 10 minutes                   |
| Packaging          | Double pouched/wrapped in CSR | Tyvek® sterilization pouch   |

Steris V-Pro Low Temperature Sterilization

| Sterilizer                                                                                                                                                                                       | Cycle                                 |
|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------|
| V-PRO® 1 Low Temperature Sterilization System                                                                                                                                                    | Standard Cycle                        |
| V-PRO® 1 Plus Low Temperature Sterilization System                                                                                                                                               | Lumen and Non Lumen Cycles            |
| V-PRO® maX Low Temperature Sterilization System                                                                                                                                                  | Lumen, Non Lumen and Flexible Cycles  |
| V-PRO® 60 Low Temperature Sterilization System                                                                                                                                                   | Lumen, Non Lumen, and Flexible Cycles |
| Devices must be cleaned, rinsed and dried in accordance with the instruction for use of the sterilizers. Refer to STERIS's Operator Manual for proper operation of the unit and cycle selection. |                                       |

Reassemble:

- 1. Snap AAV into swivel port elbow, ensuring both elbow locking tabs are fully engaged with slots of the AAV.
- 2. Press rubber cap plug onto sampling port of mask adapter.
- 3. Install swivel port assembly into face mask grooved opening. The flange of the mask adapter mates with the groove of the mask opening.
- 4. Mount headgear to face piece with the snap-in quick-release clips.


Functional Check:

- 1. Swivel port joints swivel freely; swivel port assembly is completely engaged in mask face.
- 2. AAV is installed completely in port with locking tabs fully engaged and the diaphragm flexes freely with no obstruction and functions as described in the safety features section. Do not use if it is non-functional.

Storage:

- 1. Components should be completely dry prior to storage in a sealed poly bag
- 2. Label poly bag with disinfection/clean status, mask description, date and initials.

Technical Specifications:

- 1. Operating pressure: 4 – 40 cm H2O
- 2. Mask Resistance to flow: 0.5 cm H2O at 50 L/min, 1.0 cm H2O at 100 L/min
- 3. AAV Resistance to flow (open port): 3.5 cm H2O at 50 L/min
- 4. AAV opening/closing pressure: 1.5 - 3 cm H2O
- 5. Unintentional mask leakage (leak at swivel port connections)  
Pressure (cm H2O) 10 20 30  
Leak (L/min) 5 8 10
- 6. Deadspace volume (ml)  
L(154) M(136) S(110) XS(99) P(89)
- 7. Sound pressure level: <25 dB(A) measured 15 cm from mask

Mask & Headgear Service Life

Mask, mask adapter and headgear when maintained and reprocessed per these instructions are expected to stay in service for 25 disinfection or sterilization cycles and the headgear 25 cleaning cycles. Visual inspection of the components during reprocessing will identify any defects. If defects are identified dispose of and do not use.

Operational Environment, Mask Assembly & Headgear:

- 1. Storage/Operational: Temperature: 5 – 40 °C, Humidity: 0 – 95% RH
- 2. Extreme Ambient Transport Conditions (Temporary): -40 – 70 °C, Humidity: 0 – 95% RH

Recommendations for Disposal:

Treat as conventional solid waste in accordance with local and federal regulations.

Credits:

Velcro™ is a Trademark of Velcro, U.S.A.  
Cidex™ is a Trademark of Johnson & Johnson Medical Products, Inc.  
Tyvek® is a Trademark of DuPont  
Steris® is a Trademark of Steris

**HANS RUDOLPH, Inc.**  
8325 Cole Parkway, Shawnee, KS 66227 USA  
Tel: (913) 422-7788, Toll Free: (800) 456-6695, Fax: (913) 422-3337  
Website: www.rudolphkc.com, Email: hri@rudolphkc.com

**E C Rep Ltd.**  
5 Fitzwilliam Square East, Dublin 2, D02 R744, Ireland  
Tel: +353 1 2 544 944, Email: info@ecrep.ie


**Swiss AR Services GmbH.**  
Industriestrasse 47  
CH-6300 Zug




Responsible Person in the United Kingdom  
Responsible Person Ltd  
Summit House, 4-5 Mitchell Street, Edinburgh, EH6 78D, Scotland

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Cidex™ är ett varumärke som tillhör Johnson & Johnson Medical Products, Inc.