HANS RUDOLPH, INC. 6600 SERIES V2 MASK™ ORO-NASAL CPAP/BILEVEL SINGLE PATIENT USE with VENTS & AAV INSTRUCTIONS FOR USE

Figure 1


Figure 2


Figure 3

Figure 4


Figure 5



Figure 6

Intended Use:

The 6600 V2 Mask is a disposable single-patient use, Oro-Nasal CPAP/BiLevel mask which incorporates a passive, continuous flow exhaust port. It is intended for use on adults (>30 kg) in hospitals and other clinical settings with CPAP/BiLevel machines for obstructive sleep apnea treatment, and for use with other similar ventilators that use this mask exhaust port configuration providing minimum 3 cm H₂O pressure at mask.

Caution:

Federal law restricts this device to sale by or on physician's order.

Warnings:

- Do not attempt to sterilize or disinfect mask or headgear. Sealing surfaces and headgear will lose shape under these processes causing excessive leaks. Follow cleaning instructions.
- To avoid skin irritation, never use any cleaning solutions or chemicals on this mask or headgear other than mild detergent and water.
- Inspect and clean before use, replace if not cleanable.
- Patients with facial hair may experience mask leakage.
- Do not block mask vent holes or patient circuit flow.
- Do not use with ventilators that require double-limb patient circuits.
- Not for use with CPAP/BiLevel devices which include separate patient circuit exhaust port.
- When CPAP/BiLevel device is turned off, or is not operating properly, flow through mask vent holes may be inadequate to clear all exhaled gas and some rebreathing may occur which could lead to suffocation.
- Discontinue mask use if irritation or allergic reaction develops due to mask.
- Patient failure to arouse and remove mask after vomiting could result in vomitus aspiration.
- If supplemental oxygen is used, inhaled percent oxygen will vary depending upon pressure settings, patient breathing pattern, mask size & leak rate.
- Oxygen flow must be turned off when CPAP/BiLevel system is not operating.
- Combination of this mask with devices other than for which mask is intended can alter mask safety and performance becoming unsafe or life-threatening.

Contraindications:

- Open wounds prone to infection.
- Hemodynamic or cardiorespiratory instability.
- Unconsciousness or patients unable to remove mask.
- Claustrophobia, anxiety, or other discomfort with Oro-Nasal mask.
- Facial or nasopharyngeal deformity or other inability to fit mask & seal properly.
- Excessive reflux, GI blood or other secretions.
- Impaired cough reflex, hiatal hernia, or inability to swallow or clear secretions.
- Barotrauma, upper airway obstruction or facial trauma.
- Need for ventilation or ventilator support > 12 hours per day.
- Recent facial, esophageal, or gastric surgery.
- Patients under medication with drug that may cause vomiting.
- Patients requiring immediate intubation.

Risks and Complications:

- Biosensitivity to device materials (irritation, sensitivity, or other reaction),

- Irritation/abrasion/ulceration/pain/deformity/drying (face, head, neck, eyes, mucosa, dental),
- Inhaled, aspirated, applied or handled contaminants (particulates, microbes, body fluids, chemicals),
- Infection (skin, eyes, mucosa, respiratory),
- Excessive rebreathed dead space or breathing resistance (dyspnea, atelectasis, barotrauma),
- Excessive/damaging sound level,
- Low secretion clearance (especially during upper respiratory infections),
- Inadequate or excessive mask operating pressure (barotrauma, dyspnea, atelectasis, gastric distention, flatulence),
- Adverse interactions with attached non-HRI accessories.

Safety Information:

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

SYMBOLS	
Medical device manufacturer	Needs protection from UV light sources
Do not use if the package is damaged or opened	Temperature limitations
Use-by date (first use)	Consult the instructions for use
LOT Manufacturer's batch/lot code	CE 2797 European conformity (CE) mark
REF Manufacturer's catalogue/part number	MD Indicates product is a medical device
Date of the manufacture	# Model (Series) number
Single Patient multi use	UK CA 0086 United Kingdom conformity mark
Rx Only For sale by/on physician's order only (USA)	L Size large
M Size medium	S Size small
XS Size extra small	P Size petite

Mask Components (fig 1):

- Mask face piece has plastic braces with slot openings for attachment of headgear strap mounting quick-release clips.
- Swivel port assembly comprises of sampling port, micro vent holes, anti-asphyxia valve (AAV) and detachable 22 mm (ISO 5356-1) swivel port.
- AAV (Fig. 2) snaps into swivel port assembly elbow using "locking tabs" for engagement.

Application & Fitting (fig 4):

- Determine patient mask size with optional Mask Sizing Gauge (fig 3).
- Place mask over patient's nose and mouth (1). Chin should fit into mask chin cup portion with top of mask sealing area on the bridge of the nose.
- Slide headgear over patient's head (2) with one of the lower quick release strap clips disconnected. Reconnect headgear strap clip after headgear is placed (3). Bottom straps should be positioned below ears and top straps above ears & below eyes.
- Adjust headgear Velcro™ strap tension by pulling back slightly on both the top and bottom straps (4 & 5). Adjust straps at crown of head by pulling each in opposing directions until reaching required tension and then reconnecting Velcro™ (6).
- Connect mask to patient circuit by following ventilator operating instructions.
- Turn on ventilator. If you detect leaks reposition mask or adjust strap tension. If leak continues, regardless of adjustments, try another mask size.

Safety Features Verification:

- AAV (fig 5) allows room air breathing for patient if ventilator device stops. AAV diaphragm CLOSES (>3 cm H₂O) the room air port in mask during CPAP/BiLevel ventilation. It OPENS (<3 cm H₂O) mask to room air breathing when ventilator device is off. Mask should not be used if AAV is missing, damaged or not functioning properly.

2. Mask Removal with Quick-Release Headgear (fig 6). To remove mask slide one of the lower strap quick-release clips slightly forward towards mask front. It will disconnect from mask; pull mask and headgear to opposite side or up over patient's head.

Disassembly for Cleaning:

1. Separate headgear from mask.
2. Separate swivel port assembly from mask, AAV from swivel port (squeeze two locking tabs & pull AAV out) and rubber cap plug from mask adapter sampling port.

Cleaning Mask face piece, swivel port components & headgear:

Manual Cleaning (Mask & swivel port assembly):

1. Soak components for 5 minutes in warm water with mild detergent (neutral pH).
2. Hand wash components with sponge or soft brush; then rinse with warm tap water.
3. Place components in warm water bath and agitate for two minutes.
4. Rinse in clean tap water for at least 1 minute.
5. Allow components to air dry or dry with clean, lint free cloth.

Manual Cleaning (Headgear):

1. Submerge headgear in warm soapy tap water; gently rub all areas.
2. Rinse in clean tap water for 1 minute or until all soap is removed; air dry.
3. Precaution: Do not use bleach, chlorine or alcohol based solutions to clean mask or headgear. These solutions can damage components. Direct sunlight exposure on mask and headgear components can cause deterioration and reduce product life.


Reassemble:

1. Snap AAV into swivel port elbow, ensuring both elbow locking tabs are fully engaged with slots of the AAV.
2. Press rubber cap plug onto sampling port of mask adapter.
3. Install swivel port assembly into face mask grooved opening. The flange of the mask adapter mates with the groove of the mask opening.
4. Mount headgear to face piece with the snap-in quick-release clips.

Functional Check:

1. Swivel port joints swivel freely; swivel port assembly is completely engaged in mask face.
2. AAV is installed completely in port with locking tabs fully engaged and the diaphragm flexes freely with no obstruction and functions as described in the safety features section. Do not use if it is non-functional.

Technical Specifications:

1. Operating pressure: 4 – 40 cm H₂O
2. Mask Resistance to flow: 0.5 cm H₂O at 50 L/min, 1.0 cm H₂O at 100 L/min
3. AAV Resistance to flow (open port): 3.5 cm H₂O at 50 L/min
4. AAV opening/closing pressure: 1.5 - 3 cm H₂O
5. Unintentional mask leakage (leak at swivel port connections)

Pressure (cm H ₂ O)	10	20	30
Leak (L/min)	5	8	10
6. Deadspace volume (ml)

	L(154)	M(136)	S(110)	XS(99)	P(89)
Deadspace volume (ml)					
7. Sound pressure level: <25 dB(A) measured 15 cm from mask

Mask & Headgear Service Life:

1. Mask and headgear assembly may be cleaned and reused on single-patient basis 7 days max. Do not attempt to sterilize or disinfect mask or headgear. Sealing surfaces and headgear will lose shape under these processes causing excessive leaks. Follow cleaning instructions.

Operational Environment, Mask Assembly & Headgear:

1. Storage/Operational: Temperature: 5 – 40 °C, Humidity: 0 – 95% RH
2. Extreme Ambient Transport Conditions (Temporary): -40 – 70 °C, Humidity: 0 – 95% RH

Recommendations for Disposal:

Treat as conventional solid waste in accordance with local and federal regulations.

Credits:

Velcro™ is a Trademark of Velcro, U.S.A.

HANS RUDOLPH, Inc.

8325 Cole Parkway, Shawnee, KS 66227 USA
Tel: (913) 422-7788, Toll Free: (800) 456-6695, Fax: (913) 422-3337
Website: www.rudolphkc.com, Email: hri@rudolphkc.com

E C Rep Ltd.

5 Fitzwilliam Square East, Dublin 2, D02 R744, Ireland
Tel: +353 1 2 544 944, Email: info@ecrep.ie

Swiss AR Services GmbH.

Industriestrasse 47
CH-6300 Zug

Responsible Person in the United Kingdom
Responsible Person Ltd
Summit House, 4-5 Mitchell Street, Edinburgh, EH6 78D, Scotland

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